

# Price List of Digital Hearing Aids


## w.e.f. 15 June, 2010




S.No.	MODEL	TYPE	BATTERY SIZE	MAX. BATT. LIFE (IN HRS.)	FITTING RANGE	FEATURES	MRP Rs.
1	Bravo B1	BTE	13	375	20-80 dB HL	2 Channels, 2 Bands, High Level Compression (HLC)	10,000.00
2	Bravo B11	Power BTE	13	375	40-105 dB HL		14,000.00
3	Bravo B1X	ITC	312	200	20-80 dB HL		12,000.00
4	Bravo B2	BTE	13	375	20-80 dB HL	2 Channels, 2 Bands with Wide Dynamic Range Compression (WDRC), can be programmed on 500 kHz, 1 kHz, 2 kHz, 4 kHz.	16,000.00
5	Bravo B12	Power BTE	13	375	40-105 dB HL		19,000.00
6	Bravo B32	High Power BTE	675	500	50-120 dB HL		22,000.00
7	Bravo B2X	ITC	312	200	20-80 dB HL		17,000.00
8	Bravo BXP	ITC Power	312	200	20-90 dB HL		18,000.00
9	Bravo B2-CIC	CIC	10	105	20-70 dB HL		19,000.00



10	Bravissimo BV-8	BTE	13	450	20-80 dB HL	Advance entry level Digital Hearing Aids with 3 Channels and 3 Bands with WDRC, can be programmed on 500kHz, 1kHz, 2kHz and 4kHz, WDRC, Sound Stabilizer, with adaptive features, Microphone Noise Reduction, Automatic Feedback Manager and Bravissimo open fitting excellence	27,000.00
11	Bravissimo BV-8è	BTE Open Fitting	13	450	HF<65: LF<20 dB		27,000.00
12	Bravissimo BV-18	Power BTE	13	380	40-100 dB HL		30,000.00
13	Bravissimo BV-38	High Power BTE	675	490	50-120 dB HL		31,000.00
14	Bravissimo BV-X	ITC	312	225	20-80 dB HL		29,000.00
15	Bravissimo BV-XP	ITC Power	312	200	20-85 dB HL		29,000.00
16	Bravissimo BV-CIC	CIC	10	125	20-75 dB HL		30,000.00



17	Real RE-9	BTE	13	370	10-90 dB HL	Widex Real series offers Integrated Signal Processing with 3 Channels and 3 Bands, Speech and Noise Tracer, S/N Ratio Tracer, Spatial Sound Tracer, Spatial Feedback Tracer, EDRC, Locator, Feedback Cancelling, Classic Noise Reduction, Battery Monitor, EcoTech II, Multi-Point Microphone Matching, DAI auto-detect, Telecoil, Sensogram	40,000.00
18	Real RE-9è	BTE Open Fitting	13	370	HF<75: LF<30 dB		40,000.00
19	Real RE-19	Power BTE	13	375	40-105 dB HL		42,000.00
20	Real RE-Micro	Micro BTE	10	130	10-85 dB HL		45,000.00
21	Real RE - X	ITC	312	130	10-85 dB HL		43,000.00
22	Real RE-CIC	CIC	10	150	10-80 dB HL		45,000.00



23	Flash FL - 9	BTE	13	370	0-80 dB HL	The Widex Flash series offers Integrated Signal Processing into the value sector of the market, delivering 5 Channels and 5 Bands, Speech Intensification System-SIS, Flash Locator, EDRC, Occlusion Manager, Sensogram, Feedback Cancelling, Automatic DAI detection, EcoTech II, Sound Diary and Flash open fitting excellence. Can be programmed upto 10 kHz in RIC	55,000.00
24	Flash FL - 9è	BTE Open Fitting	13	370	HF<75: LF<30 dB		55,000.00
25	Flash FL - 19	Power BTE	13	365	40-100 dB HL		57,000.00
26	Flash FL - Micro	Micro BTE	10	120	0-75 dB HL		65,000.00
27	Flash FL - X	ITC	312	195	0-80 dB HL		58,000.00
28	Flash FL - XP	ITC Power	312	195	20-80 dB HL		58,500.00
29	Flash FL - CIC	CIC	10	120	0-75 dB HL		59,000.00
30	Passion PA 105	RIC	10	125	0-80 dB HL		65,000.00



31	Aikia AK - 9	BTE	13	370	0-80 dB HL	The exceptionally smart digital hearing aids with Integrated Signal Processing, delivering 10 Channels and 10 Bands, Sound Analysis, Sound Diary with Data Logging, Spatial Feedback Tracer, Speech Intensification System, Aikia Locator, EDRC, Eco Tech II, Sensogram Multi-Point Microphone Matching, Occlusion Manager, Customized Listening Programs and Aikia open fitting excellence. Can be programmed upto 10 kHz in RIC, RC3 Remote control compatible	75,000.00
32	Aikia AK - 9è	BTE Open Fitting	13	370	HF<75: LF<30 dB		75,000.00
33	Aikia AK - 19	Power BTE	13	365	40-100 dB HL		77,000.00
34	Aikia AK - Micro	Micro BTE	10	120	0-75 dB HL		85,000.00
35	Aikia AK - X	ITC	312	195	0-80 dB HL		78,000.00
36	Aikia AK - CIC	CIC	10	120	0-75 dB HL		80,000.00
37	Passion PA 110	RIC	10	125	0-80 dB HL		85,000.00
38	Aikia RC3 Remote Control	RC3 -1 (M) RC3 -2 (B)	Lithium CR-2032	Aikia hearing aids integrates fully with the RC3 remote control. With this optional remote control you can even take "snapshots" of specific listening situations that may require special adjustments.		10,000.00 15,000.00	



39	Mind M3 - 9	BTE	13	325	0-80 dB HL	Latest technological development with 10 Channels and 10 Bands with Dual ISP, TruSound Compression, High Definition Locator, Multi-Directional Feedback Cancelling, Noise Reduction, Sensogram, Occlusion Manager, Audibility Extender, SmartSpeak, Sound Diary II, Multi-Point Microphone Matching, Service Tracker and Mind M3 open fitting excellence. CB receiver offers highest fidelity. RC4 Remote control compatible.	1,10,000.00
40	Mind M3 - 9è	Open BTE	13	325	HF<75: LF<30 dB		1,10,000.00
41	Mind M3 - 19	Power BTE	13	325	40-100 dB HL		1,10,000.00
42	Mind M3 - Micro CB	Micro BTE	10	110	0-85 dB HL		1,10,000.00
43	Mind M3 - X	ITC	312	165	0-85 dB HL		1,10,000.00
44	Mind M3 - CIC	CIC	10	125	0-80 dB HL		1,10,000.00
45	Mind RC4 Remote Control	RC4 -1 (M) RC4 -2 (B)	Lithium CR-2032	Mind M3 hearing aids integrates fully with the RC4 remote control. With this optional remote control you can even take "snapshots" of specific listening situations that may require special adjustments.			10,000.00 15,000.00



46	Inteo IN - 9	BTE	13	370	0-80 dB HL	Widex Inteo series offers Integrated Signal Processing, delivering 15 Channels, 15 Bands, High Definition Locator System, Multi-Directional Active Feedback Cancelling, Customized Listening Programmes, Speech Enhancer, Audibility Extender, Sensogram, Occlusion Manager, Sound Dairy, Noise Reduction, Multi-Point Microphone Matching and Inteo open fitting excellence. Can be programmed upto 10 kHz in RIC, RC3 Remote control compatible.	1,10,000.00
47	Inteo elan IN - 9è	BTE Open Fitting	13	370	HF<75: LF<30 dB		1,10,000.00
48	Inteo IN - 19	Power BTE	13	365	40-100 dB HL		1,10,000.00
49	Inteo IN - Micro	Micro BTE	10	120	0-75 dB HL		1,10,000.00
50	Inteo IN - X	ITC	312	195	0-80 dB HL		1,10,000.00
51	Inteo IN - XP	ITC Power	312	195	20-85 dB HL		1,10,000.00
52	Inteo IN - CIC	CIC	10	120	0-75 dB HL		1,10,000.00
53	Passion PA 115	RIC	10	125	0-80 dB HL		1,20,000.00
54	Inteo RC3 Remote Control	RC3-1 (M) RC3-2 (B)	Lithium CD-2032	Inteo hearing aids integrates fully with the RC3 remote control. With this optional remote control you can even take "snapshots" of specific listening situations that may require special adjustments.			10,000.00 15,000.00



55	Mind M4 - 9	BTE	13	370	10-85 dB HL	World's Finest Digital Hearing Aids with 15 Channels, 15 Bands Dual Chip (ISP), TruSound compression, HD Locator, Speech Enhancer, Multi-Directional Active Feedback Cancellation, Zen Music (fractal) Generator, Smart-Speak Speech Message Generator, Sensogram, Occlusion Manager, In-situ RECD Measurement Sound Diary II, Service Tracker, RC4 remote control compatible.	1,50,000.00
56	Mind M4 - 9è	Open BTE	13	370	HF<75: LF<30 dB		1,50,000.00
57	Mind M4 - 19	Power BTE	13	365	40-100 dB HL		1,50,000.00
58	Mind M4 - Micro	Micro BTE	10	120	0-80 dB HL		1,50,000.00
59	Mind M4 - X	ITC	312	195	0-85 dB HL		1,50,000.00
60	Mind M4 - CIC	CIC	10	120	0-80 dB HL		1,50,000.00
61	Passion 440	RIC	10	125	0-80 dB HL		1,60,000.00
62	Mind RC4 Remote Control	RC4 -1 (M) RC4 -2 (B)	Lithium CR-2032	Mind hearing aids fully integrate with RC4 remote control. With this optional remote control you can even take "snapshots" of specific listening situations that may require special adjustments.			10,000.00 15,000.00

B.T.E. - Behind the ear hearing aid I.T.C. - In the canal hearing aid C.I.C. - Completely in the canal hearing aid, R.I.C. - Receiver in the canal

NOTES : • All hearing aids carry a Worldwide warranty of 24 months • Prices indicated are maximum recommended retail price. • This price list supercedes all previous price lists & is subject to change without any prior notice. • Government taxes & levies will be charged extra as applicable.